






Compatibility Chart

Phone/Device	Wireless Hearing Instrument	App Required		Requires Phone Clip+	
<p>iPhone® XR, iPhone XS, iPhone XS Max, iPhone X, iPhone 8 Plus, iPhone 8, iPhone 7 Plus, iPhone 7, iPhone 6s Plus, iPhone 6s, iPhone 6 Plus, iPhone 6, iPhone 5s, iPhone 5c, iPhone 5</p> <p>iPad® Pro, iPad Air 2, iPad Air, iPad mini 4, iPad mini 3, iPad mini 2, iPad mini, iPad (6th generation), iPad (5th generation), iPad (4th generation), iPod touch (6th generation), iPod touch (5th generation)</p> <p>Apple Watch® Series 4, Apple Watch Series 3, Apple Watch Series 2, Apple Watch</p>	<p>ReSound LiNX Quattro™ ReSound LiNX 3D™ ReSound ENZO 3D™</p>	<p>ReSound Smart 3D™</p> <p>iOS 10 or later is required. Watch support requires watch OS 3 and iOS 10.</p>		<p>No</p>	
<p>Samsung Galaxy S10 Samsung Galaxy S10+ Samsung Galaxy S10e Samsung Galaxy S9 Samsung Galaxy S9+ Samsung Galaxy S8 Samsung Galaxy S8+ Samsung Galaxy S7 Samsung Galaxy S7 Edge Samsung Galaxy S6 Samsung Galaxy S6 Edge Samsung Galaxy S5 Samsung Galaxy S4</p>	<p>Samsung Galaxy Note 3 Neo Samsung Galaxy Note 4 Samsung Galaxy Note 5 DORO 8030 Google Pixel XL HTC One M8 HTC One M9 HTC 10 Huawei Mate 8 Huawei P9 LG G5 LG V20 Sony Xperia Z3</p>	<p>ReSound LiNX Quattro ReSound LiNX 3D ReSound ENZO 3D</p>	<p>ReSound Smart 3D</p> <p>Android KitKat (4.4.2-4.4.4), Lollipop (5.0 or higher) or Marshmallow (6.0.1) is required</p>		<p>- Yes for audio streaming - No for remote control functionality</p>
<p>iPhone, iPad, iPod Touch devices</p>	<p>ReSound LiNX Quattro ReSound LiNX 3D ReSound ENZO 3D</p>	<p>ReSound Relief™</p> <p>iOS 8 or later is required.</p>		<p>No</p>	
<p>Android devices</p>	<p>ReSound LiNX Quattro ReSound LiNX 3D ReSound ENZO 3D</p>	<p>ReSound Relief</p> <p>Android KitKat (4.4.2-4.4.4), Lollipop (5.0 or higher) or Marshmallow (6.0.1) is required</p>		<p>Yes for audio streaming</p>	
<p>Other Bluetooth Enabled Phones</p>	<p>Wireless Legacy ReSound</p>	<p>ReSound Control™</p>		<p>Yes for audio streaming and remote control functionality</p>	

